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
Reference Electrodes

Introduction

The XR100 and XR110 are general purpose calomel reference electrodes, fitted with a screw cap.

Preparation for Measurement

1. Remove the protection cap from the electrode and any seals covering the filling hole and the porous plug.
2. Before starting a measurement, remove the clip which closes the electrode filling hole. Remember to replace the clip at the end of measurements.
3. Check the level of the filling solution. It should be approximately 0.5 cm below the filling hole.
On delivery, both electrodes are filled with saturated KCl solution and KCl crystals.
It is always recommended to use a saturated solution of KCl. Saturation is assured when KCl crystals are visible.
4. In order to remove air bubbles trapped inside the electrode, shake the electrode holding it at its head with the porous pin down.

Maintenance

1. Electrode contamination is a major cause of faulty measurements.
2. The electrode should be rinsed with distilled water after measurements.
3. Check frequently the level of filling solution.
4. In case of deposits which cover the electrode, clean the electrode with:
 - a solution of acid (0.1M HCl, 0.1M HNO₃): mineral salt deposits, etc...
 - KS400 Pepsin in HCl Solution or RENOVO•X Xtra Strong Cleaning Solution: protein deposits (milk, cheese, serums...). Duration of treatment 1 to 2 hrs.
 - KS410 Thiourea Solution: for porous pins contaminated with sulphides or blocked by an AgCl precipitate. Duration of treatment, a few hours until the porous plug turns white.
 - RENOVO•N Normal Cleaning Solution : greasy

or oily deposits...

The porous plug of the electrode can be cleaned using a fine abrasive paper.

Storage

Between measurements: leave the electrode in KS100 Saturated KCl Solution.

Several days: seal the filling hole with paraffin film or with the electrode clip and immerse the electrode in KS100 Saturated KCl Solution using the accompanying cap.

Accessories

ΔKCL•C KCl Crystals, 15 g	S21M001
ΔKC•L Saturated KCl Solution, 100 ml	S21M002
KS100 Saturated KCl Solution, 500 ml	C20C300
○RENOVO•N	
Normal Cleaning Solution, 250 ml	S16M001
□RENOVO•X	
Xtra Strong Cleaning Solution, 250 ml	S16M002
KS400 Pepsin in HCl Solution, 250 ml	C20C370
KS410 Thiourea Solution, 250 ml	C20C380
CL111 Electrode cable with banana plug	A94L111

Specifications

Reference potential vs. standard hydrogen electrode when filled with sat. KCl							
Temperature (°C)	0°	10°	20°	25°	30°	40°	60°
mV	260	254	248	244	241	234	220

Temperature range: -10 to 60°C

Note: an electrode chain comprising pHG301/311 or XG100/200 glass electrode and a XR100/XR110 filled with saturated KCl solution will give a zero pH of approx. pH 6.65.

In case of contact with mercury, wash exposed areas immediately with plenty of water.

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RADIOMETER ANALYTICAL SAS
D31M073 • Radiometer Analytical SAS • France • 2005-12C