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
REF921
Reference Electrode

Introduction

The REF921 is a Hg/Hg₂Cl₂ (calomel) reference electrode with 1M LiCl salt-bridge solution. The electrode is recommended for measurements in non-aqueous solutions as LiCl is more soluble in organic solvents.

Preparation for Measurement

1. Remove the protection cap from the electrode and any seals covering the filling hole and the porous plug.
2. Before starting a measurement, remove the clip which closes the electrode filling hole. Remember to replace the clip at the end of measurements.
3. Check the level of the filling solution. It should be approximately 0.5 cm below the filling hole.
If necessary, refill with 1M LiCl solution. Other concentrations or other chlorides may also be used.
4. In order to remove air bubbles trapped inside the electrode, shake the electrode holding it at its head with the porous plug down.

Maintenance

1. Electrode contamination is a major cause of faulty measurements.
2. The electrode should be rinsed with distilled water after measurements.
3. Check frequently the level of filling solution.
4. In case of deposits which cover the electrode, clean the electrode with:
 - a solution of acid (0.1M HCl, 0.1M HNO₃): mineral salt deposits, etc...
 - KS400 Pepsin in HCl Solution or RENOVO•X Xtra Strong Cleaning Solution: protein deposits (milk, cheese, serums...). Duration of treatment 1 to 2 hrs.
 - KS410 Thiourea Solution: for porous plugs contaminated with sulphides or blocked by an AgCl precipitate. Duration of treatment, a few hours until the porous plug turns white.
 - RENOVO•N Normal Cleaning Solution : greasy or oily deposits...

The porous plug of the electrode can be cleaned using a fine abrasive paper.

Storage

Between measurements: leave the electrode in 1M LiCl solution.

Several days: cover the porous plug with paraffin film or with the electrode clip and immerse the electrode in 1M LiCl solution.

Long-term storage: seal the filling hole with paraffin film and place the electrode in its box and store at room temperature.

Accessories

○RENOVO•N Normal Cleaning Solution, 250 ml	S16M001
□RENOVO•X Xtra Strong Cleaning Solution, 250 ml	S16M002
KS400 Pepsin in HCl Solution, 250 ml	C20C370
KS410 Thiourea Solution, 250 ml	C20C380
CL111 Electrode cable with banana plug	A94L111

Specifications

Temperature range: 0 to 60°C

Note: the reference potential at 25°C for an electrode filled with 1M LiCl solution is approximately +40 mV with respect to a saturated calomel electrode.

In case of contact with mercury, wash exposed areas immediately with plenty of water.

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RADIOMETER ANALYTICAL SAS