Troubleshooting

First make sure that steps 1 to 4 described in the section "Checking" have been carried out correctly.

Causes	Solutions
Low sensitivity/slow response	Use fresh, quality buffers.
Calibration error	Soak the electrode for a few hours in pH 4 buffer. See "Cleaning".
pH glass membrane not hydrated enough or contaminated	Replace.
Electrode too old	Drifting potential
No stirring	Ensure reproducible stirring during calibration and measurement.
Flow of electrolyte through liquid junction inadequate	Check filling hole is uncovered. See "Preparing". Liquid junction blocked see "Cleaning".
Unstable readings	Replace the electrolyte solution. If necessary, soak the electrode tip in a 60°C water bath for 10 minutes.
Liquid junction clogged with crystals	Check connection on the meter and that the instrument is properly grounded.
Poor connection	Shield the beaker and/or add an inert salt.
Low ionic strength of sample	

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D31M105 • Printed by Radiometer Analytical SAS • France • 2004-05B

RADIOMETER ANALYTICAL SAS
72 rue d'Alsace, 69627 Villeurbanne Cedex, France
E-mail: radiometer@analytical.com Web: www.radiometer-analytical.com
Tel.: +33 (0)4 78 03 38 38 - Fax: +33 (0)4 78 68 88 12

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

GK2401B Electrode head
GK2401B-Elektrodenkopf
Tête d'électrode GK2401B

2.

Remove the paraffin film from the filling hole. After dry storage soak the electrode in water for 12 hours.

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3.


Before performing measurements, check the state of your electrode. See overleaf. For calibration, we recommend the use of our certified IUPAC pH standards.

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Checking

1.


The level of electrolyte solution (KCl•L, part no. S21M002) must be 0.5 to 1 cm below the filling hole. **Only refill with saturated KCl (KCl•L).**

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2.


At the measuring temperature, 0,5 cm of KCl crystals (KCl•C, S21M001) must be present in the solution. For measurements at temperatures >60°C, it is necessary to add extra KCl crystals.

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3.


Crystals should move freely in the solution. If necessary, shake the electrode like a thermometer or heat the electrode tip in warm water (60°C).

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4.

Check that there are no air bubbles trapped between the crystals. If necessary shake the electrode like a thermometer or heat the electrode tip in warm water (60°C).

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Storage

Between measurements: lower the electrode in distilled water.

For overnight storage: use pH 4 buffer.

For long-term storage (> 72 hours): store dry

1. Rinse the electrode with distilled water.
2. Cover the glass bulb and liquid junction with the empty protection cap. (Before reuse, soak the electrode in distilled water for approx. 12 hours).
3. Seal the filling hole with the electrode clip or paraffin film.
4. Place the electrode in its box and store in an upright position (electrode tip down).

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Disposal: Ordinary waste, no specific treatment required.

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Cleaning

Once a week, rinse your electrode with RENOVO•N solution. Unstable readings can be due to contamination, use one of the appropriate cleaning solutions. Rinse and recalibrate your electrode after each cleaning.

Protein

Use RENOVO•X (S16M002) solution for 3 minutes.

Use KS400 (C20C370) solution of pepsin in HCl for a few hours.

Grease or oil

Use RENOVO•N (S16M001) for 12 hours.

Use a solvent miscible with water (e.g. acetone or ethanol).

Sulphide or silver

Use KS410 Thiourea Solution (C20C380) for a few hours.

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